

Canine Site Requirements

Space

Aquatread will consult with customers as to the most appropriate configuration depending on the space available. Dogs will require walking access to straight ramps. Corner ramps can be supplied at additional cost for tight spaces to make access safer and more acceptable to patients.

Floor requirements

Drain in floor recommended.

Floor load bearing capacity must be adequate. A full treadmill can weigh 2000kg.

Electrical requirements

- 220/240V (or 380/400V)
- Dedicated circuit breaker of minimum 32A in the DB
- A power supply to our DB positioned at the machine, must be installed by your electrician (3 core x 4mm cable minimum, depending on distance of power supply from machine).
- As a safety feature, the system has its own earth leakage protection installed and cannot be connected to your existing earth leakage protector. Your electrician will install a separate isolator depending on municipal/council requirements.

Site access

The system is delivered to the site crated to be unpacked and installed by our technical personnel. No heavy equipment is required for installation. The machine can fit through most single doorways although double doorways are preferred.